T 024 7 669 2051

E dairy.breeding.evaluations@ahdb.org.uk



@AHDB_Dairy



Montbeliarde Production conversions

December 2024

When to use these formulas

These conversions should only be used for converting cow evaluations and bulls which do not have an Interbull evaluation. Please check with AHDB Dairy if in doubt.

The formulae provide a guide only to the expected evaluations of the relevant bull or cow in the UK.

How to convert a foreign proof using figures in the table below

Approximate UK PTA = a + (b * Foreign proof)

How to adjust the converted reliability

To take into account the loss of accuracy when converting a foreign proof into a UK equivalent, Interbull guidelines recommend that the foreign reliabilities should be adjusted by the squared correlations between the foreign and UK genetic evaluations.

When converting the foreign production proofs to UK equivalents, the foreign production reliability must therefore be multiplied by 0.8 (i.e. **UK reliability = 0.8 * Foreign reliability**)

				а		b		
SIM	milk	GBR	=	-556.93	+	0.394	*	CHE
SIM	fat	GBR	=	-17.26	+	0.336	*	CHE
SIM	pro	GBR	=	-15.88	+	0.414	*	CHE
SIM	milk	GBR	=	42.17	+	0.375	*	DEA
SIM	fat	GBR	=	1.77	+	0.349	*	DEA
SIM	pro	GBR	=	2.3	+	0.372	*	DEA
SIM	milk	GBR	=	67.77	+	0.319	*	FRM
SIM	fat	GBR	=	1.67	+	0.275	*	FRM
SIM	pro	GBR	=	3.27	+	0.301	*	FRM
SIM	milk	GBR	=	-2846	+	24.407	*	HRV
SIM	fat	GBR	=	-78.01	+	0.668	*	HRV
SIM	pro	GBR	=	-85.17	+	0.732	*	HRV
SIM	milk	GBR	=	-129.21	+	0.969	*	IRL
SIM	fat	GBR	=	-2.8	+	0.785	*	IRL
SIM	pro	GBR	=	-3.86	+	0.812	*	IRL



```
SIM milk GBR = -88.5 + 0.37 *
SIM fat GBR = -2.26 + 0.332 *
                              0.37 *
                                       ITA
                                       ITA
                                       ITA
SIM pro GBR =
                  -1.84 +
                              0.379 *
SIM \ milk \ GBR = 295.33 + 0.254 * NLD
                              0.233 * NLD
SIM fat GBR =
                   12.6 +
                              0.25 * NLD
SIM pro GBR =
                  11.68 +
SIM milk GBR = -2871.14 + 24.628 * SVN
                  -90.3 +
                             0.775 * SVN
SIM fat GBR =
                  -93.34 +
                              0.802 * SVN
SIM pro GBR =
                 218.32 +
                              0.301 * USA
SIM milk GBR =
SIM fat GBR =
                   1.68 +
                              0.319 * USA
SIM
                    1.6 +
                              0.323 * USA
   pro GBR =
```

Calculating percentage PTAs

Apply the following formula using the converted yield PTAs to obtain estimates for the percentage PTAs.

PTA Fat % = (PTA fat (kg) x 100) - (PTA milk (kg) x 3.99)(PTA milk (kg) + 6765) PTA Prot % = (PTA prot (kg) x 100) - (PTA milk (kg) x 3.43)(PTA milk (kg) + 6765)